

**ASD 220** (Version : AMOSASD220OP20231109 )

**Slide Dryer**

**Operation Manual**

## Foreword

Thank you for purchasing our ASD220 Slide Dryer. This operation manual introduces product functions, application methods and matters that need attention for safety. Please read in detail before operation to understand its performance and make full use of its functions. If you have any inquiries, please contact us, we are committed to delivering exceptional service whenever you need it.

After reading, please keep this manual properly for later reference.

## Contents

1. Safety Instruction.....	1
1.1 Safety Matters.....	1
1.2 Conditions of Installation.....	1
1.3 Delivery Inspection.....	1
2. Performance Features & Range of Application.....	2
3. Main Structure.....	2
4. Technical Parameters.....	2
5. Operation Guide.....	3
5.1 Keypad Description & Function.....	3
5.2. Instructions.....	3
5.3. Alarm Function.....	3
6. Component Layout & Electric Wiring Diagram.....	4
7. Schematic Diagram.....	5
8. Cleaning & Maintenance.....	6
9. Trouble Shooting.....	6

# 1. Safety Instruction

## 1.1 Safety Matters

Please read these brief guidelines. Failure to adhere to these instructions may disrupt the equipment's regular operation, lead to equipment damage, or pose safety risks.



Please use 220V $\pm$ 10%, 50Hz or 110V $\pm$ 10%, 60Hz AC.



Input power supply must be of good earth connection.



Keep away from any flammable or explosive materials.



Unauthorized equipment access risks high-voltage shock.



Only professional maintenance personnel may perform equipment repairs.



Use fuses with correct capacity.



Ensure the power socket and circuit can withstand at least twice the rated current.



Keep away from potential sources of interference.



Caution: Hot

- **Electrical Protection Ratings: I class, B type**
- **Shell protection grade: IP53**
- **Harmful liquid Leak-in proof degree: Normal (enclosed equipment without liquid leak-in proof)**
- **Working system: Continuous operation**

## 1.2 Conditions of Installation

- Provide at least 20cm of space around the equipment for proper heat dissipation.
- No water drop, steam, dust (including oily and airborne).
- No erosive, flammable and explosive gas and liquid.
- Stable and free from vibration.
- No interference from other electromagnetic signals.

- Operate within ambient temperature of 5°C-40°C and relative humidity under 90%.

### 1.3 Delivery Inspection

The product undergoes rigorous quality inspection before leaving the factory. However, negligence or severe impact during delivery may result in damage or partial loss, please check the followings:

1. Upon unpacking, please inspect the contents, including the equipment, operation manual, packing list, and spare parts.
2. Check the nameplate to verify that it corresponds to your order.
3. Ensure that there is no damage or loss during the delivery process.


## 2. Performance Features & Range of Application

ASD220 is a slide dryer designed for easy handling with an insertion angle of 60°. It features contactless temperature control, stable performance, a fully programmed control system, and easy, reliable operation. The equipment is designed to dry stretched and flattened sectioned tissue specimens after water bath.

### 3. Main Structure

Drying Shelf

Slide Shelf




Keypad

## 4. Technical Parameters

- Slide Capacity: 60 pieces max.
- Temperature Range: Ambient Temperature to 90°C
- Temperature Control Accuracy:  $\pm 1^{\circ}\text{C}$
- Rated Voltage: 220V $\pm$ 10% 50Hz (110V $\pm$ 10% 60Hz)
- Fuse : 220V/3A (110V/5A)
- Rated Power: 330W
- External Dimensions (mm): 324 (W) $\times$ 300 (D) $\times$ 85(H)
- Weight: 2.75kg

## 5. Operation Guide

### 5.1 Keypad Description & Function



**SET** Press this button for temperature setting; a single buzz confirms operation.


**HEAT** Press this button to start heating, red indicator lamp will light; press again to stop heating.

**▲** Adjust temperature by pressing this button; a single buzz confirms successful operation for increasing temperature.

- ▼ Adjust the temperature by pressing this button; a single buzz confirms successful operation for decreasing temperature.

## 5.2 Instructions

- Before operation, ensure every part of the equipment is in good condition.
- Plug in the equipment, switch on power at the rear, and listen for the buzzer. The LED screen will display the ambient temperature, and the green indicator light confirms normal operation.
- Press **SET** to set temperature, press ▲ ▼ to adjust temperature.
- Press **HEAT** to start heating.



Better histology . Better life from health



ASD220 Slide

D220OP20231109





## 8. Cleaning & Maintenance

- Before operation, carefully read the product operation manual and follow its instructions.
- Keep equipment clean. Clean the residuals and dirt after completing tasks.
- Do not use sharp or hard objects that may damage the worktable.
- Installation of fuse as Fig. 4.

Insert fuse (1) into installation part (2) and then place the entire installation part (2) into the socket (3).

- Ensure power is unplugged before replacing the fuse. Fuse specifications: 3A (220V), 5A (110V).




Fig. 4.


## 9. Trouble Shooting

Problems		Possible Causes	Judgment Methods	Solutions
No display		Poor connection to power supply	Ensure connection to power supply	Ensure connection to power supply again
		Power switch wasn't turned on	Turn on power switch	Turn on
		Fuse blow	Take out the fuse to check its internal resistance using a multimeter, if the value is infinite, the fuse is broken	Replace the fuse
Buzzer Alarm	E01 display	Indications of sensor damage, heating film malfunction, or solid-state relay issues		Replace the broken element
	E02 display	Digital sensor detects temperatures exceeding 100°C	Solid state relay or sensor is broken	Replace the solid-state relay or fuse with a new one
	E03 display	Digital sensor detects temperature is under -10 °C	Digital sensor is damaged	Replace the temperature sensor with a new one

# Standard Accessories List

No	Accessory Name	Qty	Notes
1	Slide Dryer	1 set	
2	Power Cord	1pc	
3	Fuse	2pcs	
4	Operation Manual	1pc	

Published by:



**AMOS SCIENTIFIC PTY.LTD**

24/632-643 Clayton Road, Clayton South, VIC 3169 Melbourne, Australia

Tel: +61-410190188

ABN 37 159 778 140

**[www.amos-scientific.com](http://www.amos-scientific.com)**

