

BETTER HISTOLOGY SOLUTION
BETTER LIFE FROM HEALTH



- Quik Processing
- Uniform heating
- Large Capacity

ADO 260 DRYING OVEN

Technical Specifications

- Inside dimension(WxDxH): 400x290x350mm
- Temperature range: 30 ~ 80°C
- Temperature accuracy : ± 1°C
- Rated voltage: 220V ± 10% 50HZ
- Rated power : 850W (Heater power: 800W)
- Dimension (W×D×H) : 500mm× 360mm× 620mm
- Tray size(W×D): 270mm × 400mm, Number: 2.
- Weight: 33Kgs

