



◆ EM400 Technical specification

Printing technology	UV Laser Marking
Print speed	2.5-4 seconds/cassette
Print resolution	2500dpi
Cassette type	Cassette with lid (one-piece) Cassette with lid assembled (two-piece) or without lid
Loading hopper	4
Max Capacity Per Loading hopper	40 one-piece cassettes/100 two-piece cassettes or cassettes without lid
Compatible system	Open to connect HIS/LIS systems
Dimension LxWxH	575mm x 228mm x 554mm
Weight	25kg
Power requirement	100-240 V AC 50/60Hz, 300VA
Other configuration	1. 10 inch touch screen; 2. Built-in bar-code scanner; 3. Bulit-in air purification system

◆ EM600 Technical specification

Printing technology	UV Laser Marking
Print speed	2.5-4 seconds/cassette
Print resolution	2500dpi
Cassette type	Cassette with lid assembled (two-piece) Cassette without lid
Loading hopper	6
Max Capacity Per Loading hopper	100 cassettes with lid assembled or without lid
Compatible system	Open to connect HIS/LIS systems
Dimension LxWxH	486mm x 313mm x 497mm
Weight	37kg
Power requirement	100-240 V AC 50/60Hz, 300VA
Other configuration	1. 10 inch touch screen; 2. Built-in bar-code scanner; 3. Bulit-in air purification system



LASER SLIDE PRINTER LASER CASSETTE PRINTER

Better histology solution Better life from health



AMOS SCIENTIFIC PTY.LTD.

Address: Factory 24, 632-642 Clayton Road,
Clayton South VIC3169, Melbourne, Australia
Email: info@amos-scientific.com

LASER SLIDE PRINTER

SM100

◆ Character

- ◆ Simple design, easy use;
- ◆ High printing speed, 3-5 seconds per slide;
- ◆ Noncontact UV laser technology;
- ◆ Excellent ability to export characters, numbers, letters, 1D & 2D codes, graphics;
- ◆ Fully compatible with LIS/HIS system;
- ◆ 11 inch LED touch screen via built-in software, also support to connect external PC;
- ◆ Built in bar-code scanner, and support to connect external scanners;
- ◆ Built in air purification system.

◆ Technical specification

Printing technology	UV Laser Marking
Print speed	3-5 seconds/slide
Print resolution	2500dpi
Slide type	Color Frosted slides
Loading hopper	1
Loading hopper max capacity	100 slides
Collection hopper max capacity	100 slides
Compatible system	Open to connect HIS/LIS systems
Dimension LxWxH	490mm x 230mm x 435mm
Weight	22kg
Power requirement	100-240 V AC 50/60Hz, 300VA
Other configuration	1. 11 inch touch screen;
	2. Built-in bar-code scanner;
	3. Built-in air purification system



LASER CASSETTE PRINTER

EM400

◆ Character

- ◆ Simple design, easy use;
- ◆ High printing speed, 2.5-4 seconds per cassette;
- ◆ Non contact UV laser technology;
- ◆ Excellent ability to export characters, numbers, letters, 1D & 2D codes, graphics;
- ◆ Fully compatible with LIS/HIS system;
- ◆ Support color setting print;
- ◆ 10 inch LED touch screen via built-in software, also support to connect external PC;
- ◆ Built in bar-code scanner, and support to connect external scanners;
- ◆ Built in air purification system;
- ◆ Compatible with various models of embedding cassettes on the markets



EM600

