



BETTER HISTOLOGY SOLUTION
BETTER LIFE FROM HEALTH



TEC 2602 SLIDE DRYER

Technical data of TEC 2602 Slide Dryer

- Table top: 324 × 272 mm
- Temperature control accuracy : ± 1°C
- Temperature range: ambient temperature ~90°C
- Capacity : can load 40 slides simultaneously
- Rated voltage: 220V ± 10% 50/60 Hz or 110V ± 10% 50/60 Hz
- Nominal power: 350W
- Outside dimensions (W × D × H): 324 × 300 × 85 mm
- Weight: about 2.75kg



Website: www.amos-scientific.com



AMOS SCIENTIFIC PTY.LTD.

Address: 7 Gabrielle Court, Bayswater North VIC 3153, Melbourne, Australia

Tel: +61-410190188 Fax: +61-3-98081503 Email: info@amos-scientific.com